


## COUNTING CHAMBER SLIDES

Designed for the quantitation of particulate liquid material, in the applications of water and urine analysis, cell and particle counting and etc. Molded from a high quality optical clear PS material, counting chamber slides are of a standardized depth, subject to strict quality control during production and before delivery, and ready for use, delivering considerable accuracy and precision and great convenience over the conventional glass counting chamber.




### Anatomy of the Counting chamber slide



- There are 10 chambers on the slide
- The area delimited by the grid is 3mm×3mm
- The grid is divided into 9 squares of 1mm×1mm
- Every square of 1mm×1mm is further divided into 9 smaller squares of 0.333mm x 0.333mm

Corresponding capacities:

- The whole delimited grid: 0.9 L
- Each of the 9 squares of 1mm x 1mm: 0.1 L
- Each of the 9 smaller squares of 0.333mm x 0.333mm: 0.0111 L



Formula:

To obtain the number of cells per  $\mu\text{l}$ :  $T_{\mu\text{l}} = \frac{n}{k \times N \times \text{CF}}$

Formula:

To obtain the number of cells per ml:  $T_{\text{ml}} = \frac{n \times 1000}{k \times N \times \text{CF}}$

Remarks:

n = total number of cells counted

k = 0.0111

N = number of small squares observed

CF

$T_{\mu\text{l}}$

$T_{\text{ml}}$

= concentration factor

= total cells present in  $1\mu\text{l}$

= total cells present in 1ml