

## XSZ-107T SERIES EDUCATIONAL MICROSCOPE

### Features:

- This model has recently undergone a style improvement, giving the instruments a very modern appearance. They are high power microscopes with a robust high quality design and finished to an excellent standard. It is an ideal instrument in biological, histological, pathological field and can be widely used in academic demonstration and research.



### Specification :

		XSZ-107T	
Viewing Head	Seidentopf Binocular Head	○	
	Seidentopf Trinocular Head	●	
Eyepiece	Extra Wide Field Eyepiece EW10×/ 22	●	
Nosepiece	Backward Quintuple Nosepiece	●	
	Backward Sextuple Nosepiece	○	
Objective	Infinite Plan Achromatic Objective	4×	●
		10×	●
		20×	○
		40×	●
		100×	●
		Objectives for Dark Field	100×
Condenser	Swing-out Condenser NA 0.9/ 0.25	●	
Stage	Double Layers Mechanical Stage 185×142/ 75×55	●	
Illumination	External Illumination, Aspherical Collector, Halogen Lamp 6V/ 30W	●	
	External Illumination, Aspherical Collector, Halogen Lamp 24V/ 100W	○	
	LED Illumination	○	
Photo Attachment		○	
Video Adapter	C Mount 1×	○	
	C Mount 0.5×	○	
Filter	Blue Filter	●	
	Green Filter	●	
Multi Viewing Attachments	Multi Viewing Attachments for 5 persons	○	
	Multi Viewing Attachments for 3 persons	○	
	Multi Viewing Attachments for 2 persons	○	
Polarization Attachment		○	
Phase Contrast Kit	Turret Phase Contrast Kit	○	
Dark Field Attachment		○	
Epi-fluorescent Attachment		○	
Temperature Controlled Stage		○	
Drawing Attachment		○	

Note: ●Standard Outfit, ○Optional