



## EUKITT® - MOUNTING MEDIUM

EUKITT® - technical advantages

EUKITT® is characterized by fast curing, crystal-clear optics with a favourable refractive index, good fluidity and low shrinkage. Therefore it is used very often in coverslipping machines.

The product is tolerant against water as well as traces of alcohol and is excellent suited for inclusion of pale colored preparations, that appear in media with lower refractive index rich in contrast, as well as those with colorless components.

## COVER GLASSES FOR MICROSCOPY



We deliver cover glasses for microscopy in all available sizes, thicknesses and in excellent quality and dimensional accuracy

- |                            |       |
|----------------------------|-------|
| 1. Cover Slip              | 22x22 |
| 2. Hemacytomer Cover Glass | 22x26 |



## H&E STAINING AGENTS

Stains are delivered in 1000 ml bottles. Each bottle is sufficient for 1500 up to 2500 specimens. They are registered as IVD and CE certified products and thus can be used for clinical diagnostics.

1. Harris' hematoxylin solution
2. 0.2% EOSIN Y solution alcoholic

**TruLaboratories**  
Corporation

14 Detroit Street, Cubao, Quezon City

+63 2 451 1964 | 451 1974

+63 2 451 1921

+63 917 8782288 | +63 925 8782288