

Agglutination Viewer



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Illuminates and provides a magnified view of the undersurface of test tube contents. Used in blood typing, cross matching, prothrombin time, Kahn and other agglutinations and serodiagnostic tests.

Brand: Trulab

Model: TRU-AV

Features:

- Easy to use side-mounted ON/OFF switch.
- Flexible neck allows for easy manipulation of light source.
- Input: 120-220V, 50/60 Hz, 18w. Output: 12VAC 500mA.

Certifications: UL listed

Reviews

There are yet no reviews for this product.