


Cryogenic Freezing Box



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Polypropylene racks hold either 0.5mL or 1.5 ml microtubes or 1.8 ml cryovials or items of similar diameter. These racks can be placed in a freezer with temperatures as low as –90°C (–130°F) and are stackable for convenient storage.

Brand: TruLab ®

Capacity:: 100 Placer

Reviews

There are yet no reviews for this product.