

Dry Bath



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

A dry bath is a type of laboratory equipment that is used to heat samples. It is often used in molecular biology, microbiology, biochemistry and genetic applications.

Brand: Dlab

Model: HB60-S

Features:

- Adjustable temperature from ambient to 60°C
- Excellent temperature control accuracy and uniformity to protect samples and ensure precision result of tests
- Timer function and continuous function are available to meet different applications
- Fast heating rate to improve the efficient workflow

- Especially suitable for sample preparation of semiautomatic biochemistry analyzer
- 5ml and 10ml mixed heating block
- Overheating protection

Brand: Digisystem

Model: DSD-100

Features:

- The more economical version between the two available models
- One block capacity
- Test tube insert into aluminum blocks which maintain temperature stability
- Easy temperature adjustment is made by turning a reference knob

Reviews

There are yet no reviews for this product.