


H. Pylori



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

H. Pylori Ab: is a immunochromatography based one-step in vitro diagnostic test for the qualitative determination of antibodies against H.pylori in human serum, plasma or whole blood.

H. Pylori Ag: is a chromatographic immunoassay for the qualitative determination of the Helicobacter pylori antigen in human stool sample.

-
- Serology Test
 - Immunochromatography based one-step in vitro diagnostic test • No instruments needed
 - Sensitivity : 100 %, Specificity : 97%

Rapid Test : H. Pylori

- Specimen : Serum, Plasma or Whole blood
- Reading Time : 10 min
- Storage: 2~30°C

- Immunochromatography based one-step in vitro diagnostic test
 - Detect Helicobacter pylori antigen in human stool
 - Specimen : Human stool
 - Sensitivity : 99.1% (106/107)
 - Specificity :98.3% (59/60)
 - No instruments needed
 - Reading time : 15 min
 - Storage : 1~30°C
- Brand:** SD Biosensor, Korea

Reviews

There are yet no reviews for this product.