

Hematocrit Centrifuge



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

This powerful unit provides high-speed separation while delivering advanced safety features. Maintenance-free, brushless motor provides quiet operation even at 14,000 rpm. Variable speed allows you to set your desired parameter and adjust during operation.

Brand: DragonLab and Scilogex

Model: DM1424

Capacity: 24 Placer

Features:

- Speed range of 200-14000rpm
- With hematocrit rotor and 2mlx24 rotor
- User-friendly large LCD display provides all visual information
- Quiet and stable operation
- Brushless DC motor drives quickly

- Dual door interlock design, more safe

Reviews

There are yet no reviews for this product.