

Laboratory Thermometers



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Laboratory thermometers are used to measure temperatures or temperature changes with a high degree of precision.

Brand: TruLab ® and Assistant (Germany)

Types:

- Alcohol Filled
 - With red indicating liquid and white back for readability. Length 305 mm.
 - Celsius Range: -10 - 110°C
- Ref Thermometer
 - Hangs or stands.
 - Sturdy 2 inch diameter stainless steel construction.

- Celsius Range: -30 - 30°C
- Digital Thermometer
 - Indoor & outdoor temperature and humidity display in C° or F°
 - Memory to store/recall the max/min temperature and humidity
 - Clock and date display function
 - Temperature range: -10°C to 50°C
 - Humidity range: 20% ~ 99%RH
- Room Thermometer
 - Support is made of Beech Wood
 - Length: 20cm
 - Celsius Range: -30 - 50°C

Reviews

There are yet no reviews for this product.