

Malaria P.f/PAN



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Malaria P.f/PAN: An immunochromatography based one-step in vitro diagnostic test for the qualitative determination of P.falciparum specific HRP-II (Histidin Rich protein II) and Plasmodium specific pLDH (Plasmodium lactate dehydrogenase) in human whole blood.

Brand: SD Biosensor, Korea

General Information:

- Malaria P.falciparum and the other plasmodium with separate test line at one time
- Differential detection of P.falciparum and the other plasmodium
- Sensitivity : 99.6% for P.f, 95.9% for Non-P.f
- Specificity : 99.3%
- No instruments needed
- Specimen : whole blood
- Reading Time : 20 min

- Storage : 1~30°C
- WHO PQ Approved

Reviews

There are yet no reviews for this product.