

## Multipurpose Centrifuge



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

The Multipurpose Clinical Centrifuge is ideal for the centrifugation of samples in clinical chemistry, cytology, research laboratories, industrial laboratories, and the like. This centrifuge is suitable for tubes with caps, blood collection or urine tubes. Swing-out bucket rotors and adapters can be compatible with standard tubes up to 100ml at low speed.

**Brand:** Dlab

**Model:** DM0636

### Features:

- Speed can be set from 300-6000 rpm
- Maximum capacity of the centrifuge: 100ml x 4
- 9 Acceleration / 10 Brake speed settings to ensure optimal separation
- A wide range of accessories (fixed angel and swing-out bucket rotors)
- Automatic rotor identification
- Timers can be set up to 99 minutes for short run or continuous run
- 9 procedures stored in memory makes the operation much simplified
- Maintenance-free brushless DC motors

### Reviews

There are yet no reviews for this product.