

## Parafilm M



Rating: Not Rated Yet

### Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

A unique thermoplastic, self-sealing, moldable and flexible film for numerous uses in the typical laboratory. It holds moisture loss to a minimum. It offers excellent protection for the contents of tubes, flasks, and petri dishes.

**Temperature Limit:** -45--50 °C temp. range

**Size:** 4 inches x 125 feet

### Reviews

There are yet no reviews for this product.