

Rotator



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Rotators are widely used in VDRL, RPR, and EIA lab test procedures.

Brand: Digisystem

Type: Digital (DSR-2100D) and Analog (DSR-2100D)

Features:

- Stable and quiet
- With variable speed, timer control and continuous operation

Features Digital Rotator:

- Digital display shows the present speed and remaining time
- Timer can be set from 0 to 99 hours and 59 seconds

Reviews

There are yet no reviews for this product.