

## Shaker



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

The Orbital and Rocking Shakers provide a smooth continuous motion for uniform mixing. Display of time and speed can be viewed through the double LCDs. A wide range of platform is available for use with a variety of vessels.

**Brand:** Dlab

**Motion:** Orbital

**Model:** SK-0180-S

### Features:

- Maximum load capacity is 3kg with platform
- Speed range of 40-200 rpm
- Provides continuous or timed operation

**Motion:** Rocking

**Model:** SK-R1807-S

**Features:**

- Maximum load capacity is 10kg with platform
- Speed range of 10-70 rpm
- Rocking motion angle 9°

**Reviews**

There are yet no reviews for this product.