


### Slide Staining System



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

#### Description

The Slide Staining System is recommended for manual hematology slide staining. The structural base is made of high quality stainless steel material, resistant to commonly used solvents while the reservoirs are made of POM material resistant to high temperature. The reservoirs are compatible with all standard microscope slides and can be used in the microwave oven.

**Capacity:** 3 and 12 Reservoirs

#### Reviews

There are yet no reviews for this product.