

## Specimen Seal up Racks



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

Seal Up Rack: A safe solution for blood specimen post processing.

**Features:**

**Isolated Seal**

» Avoids cross infection, ensuring safety of specimen, personnel and environment.

**Multifunctional**

» It can be separated into two parts. The bottom part can be used as universal tube rack in laboratory. No need to transfer repeatedly.

**Number Marks**

» Identify specimen easily with numbers (1~100) marks on the rack.

**One Step Operation**

» Seal up specimen in one step, improving work efficiency.

### Stackable

» Positioning groove and ridge design allows stable mounting of racks on top of each other.

### Reviews

There are yet no reviews for this product.