
Syphilis



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

Ask a question about this product

Description

A qualitative screening test kit for the detection of antibodies to Treponema pallidum in human blood.

Brand: SD Biosensor, Korea

General Information:

- Treponema pallidum Antibody Test
- Immunochromatography based one-step in vitro diagnostic test
- No instruments needed
- Sensitivity: 100% (169/170)
- Specificity: 100% (199/200)

Specimen : Serum, Plasma or Whole blood
Reading Time : 5~20 min
Storage : 1~30°C

Reviews

There are yet no reviews for this product.