

## Urine and Stool Collection



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description


TruLab®oratories offers an extensive range of urine and stool containers. The containers are made entirely of polypropylene and polysterene, are chemical-resistant even to hydrochloric acid, and are metal-, latex-, zinc-, and fluorescence-free. Containers can be subjected to freezing, thawing, and EtO gas sterilization without causing changes in materials or physical appearance.

**Brand:** TruLab®

**Accessory:** Specimen Container Rack

### Types:

- Sterile
- Non Sterile
- Volume Capacity: 30, 60, 90, and 150 ml



This 24 placer, stainless steel rack keeps specimen containers organized and in place, minimizing risks of spilling and contamination.

### Reviews

There are yet no reviews for this product.