Vortex





Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

Ask a question about this product

Description

The Vortex is a device used in laboratories to mix liquid components in tubes and resuspend cells. It is widely used in biological and chemical analysis applications.

Brand: Dlab

Model: MX-S

Features:

- Touch operation or continuous mode
- Variable speed control from 0-2500 rpm
- Used for various mixing applications with optional adapters
- Specially designed vacuum suction feet for body stability
- Robust aluminum-cast construction

Reviews

There are yet no reviews for this product.