

LIQUID URINE CONTROL LEVEL 1 AND 2Simply validate visual and urinalysis with Insight Liquid Urine Controls.



RELIABLE

- Use with Insight® Urinalysis Reagent Strips and Urine Analyzer for optimum quality control
- Validate urinalysis results and prevent procedure errors results for ASC

QUICK AND CONVENIENT TESTING

• Easily drop the control solution onto each reagent pad using the dropper tip bottle.

PACKAGING AND SHELF LIFE

- Ready-to-use without dissolving in distilled water
- 18 months shelf life for unopened controls at 2-8 °C
- Dropper Tip Bottles -available in separate positive and negative levels
- Dropper tip bottles provide efficient use of the control solution
- Control can be used up to 40 times within 20 days at room temperature

SPECIFICATIONS	
Product Name	Liquid Urine Control
Test Parameters	LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU, ASC, ALB, CRE, CA (13)
Solution Detection Level, Level 1	Negative: LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU, ASC, ALB, CRE, CA
Solution Detection Level, Level 2	Positive: LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU, ALB, CREA, CA; and Negative ASC
Compatible Urine Strips	Insight® Urinalysis Reagent Strips
Reading Time / Stability	Refer to insert
Storage Temperature	2-8℃
Opened Control Stability	24 Months
Memory	30 days at 15-30°C or until the expiration date at 2-8°C
Maximum Test Compatibility	20 or 40 tests/bottle



