



MULTIPRO SYSTEM

FAST AND SIMPLE PROCEDURE

- One step operation (20 seconds hands-on procedure)
- Fully automatic assay by analyzer, total assay time ≤ 6 minutes
- Embeds the code into every cartridge, presented as barcode
- 7 inches large touch screen
- On board printer

ACCURATE AND RELIABLE RESULT

- NGSP certified, IFCC certificate on going
- No interference from Hb variants, like HbE, HbC, HbS...
- Immunoturbidimetry Technology

ROOM TEMPERATURE STORAGE OF CARTRIDGE OPTIONAL

- 2-8°C 36-46°F: 24 months Shelf life
- 2-30°C (36-86°F); 3 months anytime before expiry date
- The cartridge has integrated all necessary reagents and sampling device

FEATURES	HEMOGLOBIN A1C TEST (HBA1c)
Methodology	Immunoturbidimetry
Test Time	≤ 6 minutes / test
Specimen	1 μ L Capillary or Venous Blood
Testing Range	3.8-15.0% NGSP, or 17.8-141.0 mmol/mol IFCC
Memory	5,000 records for patient specimen and QC test
Application	Only use On Call® Multipro Test kits and Controls
Printer	Inner Thermal printer; external printer interface
Data Transport	RS232, USB, Internet
Language	English (default in US and select countries) and additional language (s)
Operating Conditions	15-30°C (59-86°F); 10-90% RH
Power Source	240 VAC
Power Consumption	60W
Screen	TFT 1024 x 600 touch screen
Information Identification	Scans QR code automatically
Weight	2.3 kg. without power supply and printing paper
Dimensions	257mm x 313mm x 158mm

OTHERS:

HBA1c Liquid Control

- Level 1: 5 x 1 μ L
- Level 2: 5 x 1 μ L
- 10 pcs. per box

HBA1c Test Cartridge

- 10 tests per kit