

TABLE TOP CENTRIFUGE

Table Top Centrifuges are widely used in many fields such as biology, medicine, pharmacology, laboratory, and other regions. Table Top Centrifuges are small, compact, strong and easy and safe for operation. The max. Speed of Table Top Centrifuge is 3000 rpm at 60 Hz.



FEATURES:

- The design is up-to-date, beautiful outlook and durable.
- The speed and time of Table Top Centrifuge are adjustable.
- ABS outer shell can resist collision.
- Brushless, free maintenance.
- Table Top Centrifuge could accommodate 2ml~15ml centrifuge tubes without changing the rotor.
- With safety switch, when cover is opened, the motor can be stopped.
- Speed is linear and stable.
- Auto-balancing system.
- Quality certificate: ISO 9001, 13485

ROTORS



TUBE ADAPTER



SPECIFICATION

MODEL	Max. Speed (RPM)	Timer (Min.)	Speed & R.C.F. Display	Manual Brake	Motor	Weight (kg.) (N.G/ G.W.)	Dimension (mm) (W x D x H)
DSC-200T	3,000 at 60Hz 2,500 at 50Hz	0~30 min. mechanical	None	None	Brushless Induction	3.5 / 4.5	280x280x255
Power Supply AC 220V 50/60Hz available							

SPECIFICATION FOR ROTORS

Model	Speed (RPM)	Type	Max. RCF (g)	Capacity	Standard Accessory
AR-1506	300~3,000	45° angle	1,066	6x (10~15) ml	(10~15) ml tube adapter
AR-1508	300~3,000	45° angle	1,066	8x (10~15) ml	(10~15) ml tube adapter
AR-1512	200~3,400	45° angle	1,680	12x (10~15) ml	(10~15) ml tube adapter
AR-1524	200~3,400	45° angle	1,680	24x (10~15) ml	(10~15) ml tube adapter