



## DRY BATH

A dry bath is a type of laboratory equipment that is used to heat samples. it is often used in molecular biology, microbiology, biochemistry and genetic applications.

### FEATURES:

- A wide range of temperature control up to 120°C
- Timer function support is available
- Various types of optional blocks are available
- Overheating protection
- Lid for protection and heat preservation



Model: HB-120S

SPECIFICATIONS	HB-120S
Function	Heating
Temperature Range	Room Temperature - 120°C
Temp. Control Accuracy [at 37°C]	±0.5°C
Temperature Uniformity	±0.5°C
Max. Heating Rate	5.5°C/min.
Timer	1 min. - 99h 59 mins.
Screen	LED
Overheating Protection	140°C
Adapter Block Material	Aluminum
Voltage, Frequency	220-240V, 50Hz/60Hz
Power	160W
Dimension {D x W x H}	175 x 290 x 85mm
Weight	3 kg.



**TruLaboratories**  
Corporation

+632 82-TRULAB (878522)  
+632 8287.3683  
+63917.878.2288 • +63925.878.2288  
[www.trulaboratories.com](http://www.trulaboratories.com)