

LAB INCUBATORS DSI-300D

Lab incubators are necessary equipment for any laboratory conducting cell culture and tissue culture work. Incubators protect cells from changes in temperature, humidity, carbon dioxide and oxygen.

FEATURES:

- Simple keypad input allows easy temperature setting.
- LED digital display enables users to monitor the chamber temperature at any given moment.
- Visual alarm indicator alerts users of abnormal conditions if the chamber temperature exceeds the setting points by 10°C.
- Time setting can be continuous mode or timed mode (999 hours or 999 minutes).
- The temperature can be controlled and maintained to 70°C.
- Shelves can be adjustable.
- Double doors allow sample viewing from the acrylic door.
- Quality certification: ISO 9001, ISO 13485



SPECIFICATIONS	DSI-300D
Convection	Natural convection
Working Temperature	Ambient +5°C~70°C
Time mode	Continuous mode / timed mode (1~999 minutes)
Capacity (Liters)	30 Liters
Chamber Dimension (mm) (W x D x H)	325 x 310 x 315
Power Watts	300W
Overall Dimension (mm) (W x D x H)	425 x 420 x 610
Accessory	2 Shelves
Power Supply	220V 50/60Hz available