



MULTI PURPOSE CLINICAL CENTRIFUGE DM0636

The Multipurpose Clinical Centrifuge is ideal for the centrifugation of samples in clinical chemistry, cytology, research laboratories, industrial laboratories, and the like. This centrifuge is suitable for tubes with caps, blood collection or urine tubes. Swing-out bucket rotors and adapters can be compatible with standard tubes up to 100ml at low speed.

FEATURES:

- Speed can be set from 300- 6000rpm
- Maximum capacity of the centrifuge: 100 mL x 4
- 9 Acceleration /10 Brake speed settings to ensure optimal separation
- 9 procedures stored in memory makes the operation much simplified
- A wide range of accessories (fixed angle and swing-out bucket rotors)
- Automatic rotor identification
- Timers can be set up to 99 min. for short run /continuous run
- Maintenance-free brushless DC motors



Model: DM0636

Capacity: Rotor 24 Placer & 36 placer

SPECIFICATIONS	DM0636
Speed Range	300-600 rpm Increment: 10rpm
Max. RCF	4300 x g
Speed Accuracy	±20rpm
Max. Capacity	100mLx4
Run Timer	30sec-99mins / continuous
Motor	Brushless DC
Display	LCD
Program Memory	9
Safety Devices	Door Imbalance protection Over speed & Over temperature detection
Acceleration / Deceleration Time	9↑ / 10 ↓
Power	660W; Single-phase AC220V-240V, 50Hz/60 Hz, 6.3A
Dimension [W x D x H]	445 x 579 x 269 mm
Weight	36 kg.
Certification	CE
Advance Features	Speed / RCF switch; Rotor identification