


DLAB

MICROPETTE MULTI CHANNEL



8 Channel Pipette



12 Channel Pipette

Comfortable finger support

Push button for volume adjustment

Low force tip ejector for easy operation

Spring loaded tip cones for easy cleaning and maintenance

High chemical-resistant materials to secure long lifetime

FEATURES:

- 8 and 12 channel pipettes are appropriate for 96 well plates.
- Dispensing head rotates for effortless pipetting convenience.
- Individual piston and tip cone assemblies allowing easy repair and maintenance.
- Compound material - made tip cone secures high sealing performance.
- Compatible with most universal tip brands.

121°C
HALF AUTOCLAVE

MICROPETTE (Mechanical)



FEATURES:

- Pipettes cover a volume range from 0.1 µL to 10mL
- Ergonomic design provides excellent operating experience
- Large display window allows for easy volume identification
- Easy calibration and maintenance
- Each MicroPette supplied with an individual calibration certificate according to ISO8655

MICROPETTE PLUS (Autoclavable)



FEATURES:

- Fully autoclavable
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- The pipettes cover the volume range from 0.1 µL to 5mL
- Easy calibration and maintenance
- Manufactured from innovative materials
- Each MicroPette Plus supplied with an individual calibration certificate according to ISO8655

TYPES	VOLUME
Single Channel Adjustable Volume	0.1 µl-2.5 µl, 0.5 µl-10 µl, 2 µl-20 µl, 5 µl-50 µl, 20 µl-200 µl, 50 µl-200 µl, 100 µl-1000 µl, 1000 µl-5000 µl
Single Channel Fixed Volume	5 µl, 10 µl, 20 µl, 50 µl, 100 µl, 500 µl, 1000 µl
Auto Clavable	10 µl-100 µl, 100 µl - 1000 µl

PIPETTOR STAND

The linear and round stands have been designed to fit perfectly our pipette range, MicroPette and MicroPette Plus Autoclavable Pipette. The stands are convenient to hold up to 6 pipettes, both single and multi-channel.



ROUND

- Rotates smoothly



LINEAR

- Static



TruLaboratories
Corporation

+632 82-TRULAB (878522)
+632 8287.3683
+63917.878.2288 • +63925.878.2288
www.trulaboratories.com