



TUBE ROLLER

The DLab Tube Rollers provide a gentle and highly efficient rolling motion, ideal for mixing blood samples, viscous substances, and liquid-solid suspensions. Used in variety of applications such as prevention of blood coagulation, immune precipitation and many more.

FEATURES:

- Speed range of 0-70rpm
- Six roller design for gentle rolling and easy to clean



SPECIFICATIONS	MX-T6-S
Amplitude	24 mm
Motor Type	DC Motor
Rocking Motion	Rocking and rolling
Amplitude	24 mm
Max Load Capacity	4 kg
Number of Rollers	6
Roller Size [Length]	280mm
Speed Range	0-70 rpm
Speed Display	Scale
Operation Mode	Continuous
Voltage	220V-240V, 50/60 Hz
Power	25W
Dimension [D x W x H]	450 mm x 260 mm x 120 mm
Weight	4.5 kg.
Permissible Ambient Temperature & Humidity	5-40°C, 80% RH
Protection Class	IP21