

WATER BATH

A water bath is a device that maintains water at a constant temperature.



Model: DSB-500D
Type: Circulating
Capacity: 10L



Model: WB-500D
Type: Unstirred
Capacity: 10L

FEATURES:

- DSB-500D are circulating water bath which have immersion pump
- WB-500D are unstirred water bath
- Large bath opening let you immerse samples conveniently
- Inside stainless steel tank
- Long life immersion stainless pipe heater
- Quality certificate: ISO 9001, ISO 13485

SPECIFICATIONS:

Model	DSB-500D	WB-500D
Stirring Type	Circulating	Unstirred
Capacity	10 Liters	10 Liters
Inner Tank Dimension (WxDxH)	240x300x150 mm	240x300x150 mm
Temperature Setting / Read out	Digital	Digital
Temperature Control	PID	PID
Temperature Sensor	PT 100 Ω	PT 100 Ω
Stability	$\pm 0.2^{\circ}\text{C}$	$\pm 0.2^{\circ}\text{C}$
Working Temperature Range	Ambient + 5~100 $^{\circ}\text{C}$	Ambient + 5~100 $^{\circ}\text{C}$
Power Watt	1200W	650w
Optional Accesory	Flat Lid, Gable Cover, Test Tube Rack	Flat Lid, Gable Cover, Test Tube Rack
Dimension (mm) (W x D x H)	275 x 335 x 190	290 x 400 x 240
Weight (kg) (N.W./G.W.)	8.0 / 9.5	7.0 / 8.5
Power Supply	AC 20V 50/60Hz available)	AC 20V 50/60Hz available)