

PIPETTE SHAKER PS-600T

The Pipette Shaker PS-600T can easily achieve the uniform mixing of the chemicals with blood while they are in the diluting pipettes.

FEATURES:

- Has a pipette holder made of stainless steel spring clips.
- Provides continuous or timing operation (0-15 minutes)
- Can hold up to 6 pipettes
- Quality Certification: ISO 9001, ISO 13485



Model	Capacity	Switch	Fixed Speed	Weight (KG) (N.W. / G.W.)	OVER ALL DIMENSIONS (WxDxH)
PS-600T	6 Pipettes	3 Positions (ON/OFF/0-15 min. timer)	2,500 rpm	1.8/2.2	150 x135x160mm

Power Supply: AC (220V 50/60Hz available)