

ALCOHOL LAMP



An Alcohol Burner is a piece of laboratory glasswares used to produce an open flame. Available in 30ml, 60ml, 100ml and 150ml.

CARBOY WITH SPIGOT



Carboy is a leak proof polypropylene screw closure molded in spigot which provides strength and leak proof operation. Available in 5 liters and 10 liters.

BEAKER WITH AND WITHOUT HANDLE



Glass (Without Handle)

Plastic (With Handle)

A beaker is a cylindrical container used to store, mix and heat liquids in laboratories. Beakers usually have a flat bottom and a lip around the top. Available in 50ml, 100ml, 150ml, 250ml, 500ml, 1000ml, 2000ml for Glass. 250ml and 500 ml for PP material.

BOILING FLASK



Boiling flasks are used for heating and boiling of liquids. The flasks are designed to have round bodies with long, thin necks to facilitate uniform heat distribution and minimize evaporation. Its typically made of borosilicate glass, which is resistant to thermal shock and chemicals. Available in 150ml, 250ml, 500ml and 1000ml.