

Clean safely glasses and googles, microscope lens and more with no scratching.



APPLICATOR STICK

Disposable wooden applicator sticks with uniform cylindrical shape and smooth finish. They can be used for routine spot test sampling and drop transfers or wherever handy, disposable mixing is required.

Quantity: 1,000 pcs. per box

Type: Non-Sterile



WASH BOTTLES AND RACK

Wash bottles saves time, helps avoid mistakes, and increase safety. Bottles are durable, chemical-resistant, squeezable and easy to clean. The leak proof wash cap assembly will deliver a single drop or a needle thin stream, with finger-tip squeeze control flow. Available in 250 and 500 ml.

Wash bottle rack transports wash bottle. Wash Bottle rack offers safe, convenient, lightweight transportation that is sturdy and lightweight.



PH Paper is used to determine if a solution is acidic, basic or neutral. This is determine by dipping part of the paper into a solution of interest and watching color change.



PARAFILM M

Parafilm M all-purpose laboratory film covers and shields your productsfrom moisture while allowing gas permeability. Keep moisture loss to a minimum and protect the contents of your laboratory items with a flexible covering for test tubes, beakers, petri dishes and more. Convenient, flexible and self-sealing Parafilm M laboratory film is an essential item in health care, pharmaceutical and research labs.



+632 82-TRULAB (878522) +632 8287.3683
+63917.878.2288 •+63925.878.2288
www.trulaboratories.com