





WBC PIPETTE

This is a bulb pipette having a long stem with a capillary bore and a pointed tip. The bulb contains a white bead inside. A small rubber tube provided with a mouth piece is connected to the small narrow portion, above the bulb for sucking blood and fluid into the pipette.

RBC PIPETTE

This is a bulb pipette having a long stem RBC pipette is graduated to give dilution 1:00 or 1:200. It has bottom marked with 0.5 and 1 and the top is marked with 101. It has a round bulb containing red bead. A rubber tube is attached to the top for sucking.



SAHLI PIPETTE

Sahli pippetes usually used as Hemoglobin pipettes for dosing the blood with capacity of 20 µl.



DISPOSABLE TRANSFER PIPETTE

- Made from flexible polyethylene.
- 100% inert for accurate sampling
- Transparent wall permits easy viewing of pipette contents
- Clearly visible with either clear or turbid samples.
- Available in 1 ml and 3 ml Non-sterile.
- Also available in 3 ml Sterile.



This is a bulb pipette having a long stem Flexible, variable-capacity bulbs that are squeezed to create a vacuum or apply pressure to the liquid contents of pipets to collect transfer and dispense liquid.



WARDING TRAY

Blood Collection Tray is light weight and easy to clean. Deep tray design helps to reduce spills. Also available Warding tray with tube slots.



PHLEBOTOMY ARM SUPPORT

- Standard Coating option has protective vinyl layer that wipes clean and won't absorb fluids
- Antimicrobial Coating option has extra-protective vinyl layer with a silver ion antimicrobial agent built in; wipes clean and won't absorb fluids
- Keeps the elbow straight and elevated
- Provides comfort and relaxation for the patient

