

## **MICROHAEMATOCRIT TUBES**

Vitrex® micro-hematocrit tubes are quality product that conform to international standard 12772. The tubes are supplied fire polished at the colour-coded end, which we recommend should be sealed with Vitrex Sigillum Wax. Fire polishing reduces wear on centrifuge seals. Supplied with colour-coding on the packaging and on the tubes in plain (blue colour), or with uniform sodium heparanisation (red colour).

## **AVAILABLE IN:**



## **ALCOHOL SWABS**

Vitrex<sup>®</sup> Alcohol Swabs are an easy and convenient way to cleanse the skin and help fights bacteria at the site of application. Use as antiseptic skin cleanser before injections and blood sampling.

Each alcohol swabs is individually sealed and disposable for convenience and safety.

SIZE	PACKING	CONTENT
30 x 60 mm	100 pcs per box	70% Isopropyl Alcohol

