## TruScreen

## **PREGNANCY TEST CONTROLS**

TruScreen Human Chorionic Gonadotropin (hCG) Positive and Negative controls are prepared from human based urine. These hCG controls can be used to monitor the performance of qualitative immunochromatic hCG detection devices.

- TruScreen hCG Liquid Controls Urine, Positive, 5 ml vial
- TruScreen hCG Liquid Controls Urine, Negative 5ml vial



## **DRUG TEST LIQUID CONTROLS**

TruScreen Drug Test Liquid Control is a In Vitro Diagnostic IVD device, that is intended for use as quality control urine to monitor the precision of laboratory urine toxicology testing procedures for the analytes listed below:

TEST	SAMSHA MANDATED	NEGATIVE	2X
THC	Delta-9-THC-COOH	0	100
COC	Benzoylecgonine	0	600
PCP	Phencyclidine	0	50
OPI/MOR	Morphine (Low Opiate) Morphine (High Opiate)	0 0	600 4000
AMP	d-Amphetamine	0	2000
MET	d-Methamphetamine	0	2000

TEST	SAMSHA MANDATED	NEGATIVE	2X
BARB	Secobarbital	0	600
BZO	Oxazepam	0	600
MTD	Methadone	0	600
MQL	Methaqualone	0	600
PPX	Propoxyphene	0	600
TCA	Notriptyline	0	2000
MDMA	MDMA*	0	1000
OXY	Oxycodone	0	200
BUP	Buprenophine	0	20





- TruScreen with MDMA, OXY & BUP, 2x Cutoff, 5ml vial
- TruScreen Negative Control Urine, 5ml vial

