

LANCETS



Lancets are small needles that are used to prick the skin to obtain a blood sample. Available in Twist and Flat types.

#### LANCING DEVICE



A device that holds a lancet firmly, and when triggered, moves the lancet linearly ahead to prick the skin in a controlled manner. The purpose is to obtain a small amount of capillary blood for testing.

### SAFETY BLOOD LANCETS

- The needle retracts automatically and permanently into the safety lancet after use, avoiding the risk of sharp injuries.
- Safety lancets are designed for single use, avoiding the risk of cross-infection between patients.
- Each safety lancet is designed to puncture & penetrate to predetermine depth that is safe and specific to patient type.
- Safety lancets offer clear safety advatages in a clinical environment for patients, clinicians and cleaning personnel.
- Available in Gauge 21 and 26

# FLASHBACK NEEDLES



Visible Flashback Blood Collection Needle is made of stainless steel, needle point is tri-bevel cut which minimize the pain while puncture into vein, connect body at middle is transparent, which make vein blood flow visible after puncture. Available in gauge 21 and 22.

# **SAFETY BLOOD COLLECTION SET**



Safety blood collection sets are designed to increase the safety of routine venipuncture procedures. The winged blood collection needle has a safety mechanism that is activated after blood collection as the needle is pulled out of the vein-for maximum protection against needlestick injuries. Available in Gauge 21, 23 and 25.

# MULTISAMPLE NEEDLES

Multisample Blood Collection Needles permit several samples to be taken with a single puncture. Sharp, smooth edges help reduce pain upon penetration, and provide easy connection to rubber stoppers. Available in gauge 21, 22 and 23.



### **NEEDLE ADAPTER**



These needle adaptors allow for economical extension Single and multiple-use devices used in conjunction with blood collection tubes and needles to draw and transfer patient blood samples in healthcare settings.

## **ADJUSTABLE TORNIQUET**

A tourniquet is a device which applies pressure to a limb or extremity in order to limit-but not stop the flow of blood. It may be used in emergencies, in surgery, or in post operative rehabilitation.



Available Colors: Light Blue, Dark Blue, Plum, Red & Green



+632 82-TRULAB (878522)

+632 8287.3683 +63917.878.2288 • +63925.878.2288 www.trulaboratories.com