

SLIDE TRAY



Storing slides horizontally lets you check your slides visually. Colors available in Blue, pink, and white. With & without transparent cover.

SLIDE BOX



Hinged ABS plastic box has metal clasp, index in lid and cork liding. Number slots holds 25-100 slides. Stackable.

TRAY DISPENSER



Easy-to-use top-loading dispenser accepts standard $25 \times 75 \text{ mm} (1'' \times 3'') \text{ slides}.$ Capacity: 91 Slides.

SLIDE MAILER







Double Capacity



Triple Capacity

Economic & resealable, made of extra-strong, clear polypropylene, these slides mailers keep slides secured during transport.

SHARPS CONTAINER

Sharps container systems provide solutions for the safe and effective containment and disposal of contaminated needles, syringes and other sharps. They are stackable, puncture-resistant and autoclave-and incinerator-safe. They are featured transparent, locking lids and needle removal ports. The tortous path systems restrict hand entry.



0.7 Liter



1 Liter



5 Liters



15 Liters



23 Liters



Laboratory Waste Can 20 Liters capacity



Needle Destroyer with Cutter





www.trulaboratories.com