


Virus Collection and Transport Kits

- The complete solution designed for the collection, preservation, and transportation of viral specimens.
- Ready to use and easy to tear package, effectively avoiding cross contamination.
- Supplied with biohazard specimen bags to ensure safe transport of samples.

6 EASY STEPS:



Ordering Information	
Kit Components	Quantity
Nasopharyngeal Swab, Sterile (Flocked)	1 piece
Oral Swab, Sterile (Viscose)	1 piece
10ml tube with 3ml VTM Medium Hank's Balanced Salt Solution, Bovine Serum Albumin, Glucose HEPES Buffer, Penicillin-Streptomycin, Amphotericin B, Phenol Red	1 piece
Biohazard Bag	1 piece





- ¶ 14 Egea St. Cubao, Quezon City
 8451 1964 | 8451 1974 8451
- **2** 8451.1964 | 8451.1974 **3** 8451.1921 **3** 0917 878.2288 | 0925 878 2288
- www.trulaboratories.com