



DENGUE TEST KITS

Wondfo Combined Diagnostic Kit for Dengue NS1 Antigen and IgG/IgM Antibody is a rapid test for the detection of dengue NS1 antigen and IgG/IgM antibodies in human whole blood, serum or plasma as an aid in the diagnosis of dengue infection. For in vitro diagnostic use only. For professional use only.

GENERAL INFORMATION:

- Accurate
- Sample: Serum/Plasma/ Whole Blood
- Early detection of Dengue Virus
- Distinguish primary and secondary infection
- Visual result, simple and safe to use




EXCELLENT CLINICAL COMPARISON RESULT


ITEMS	SENSITIVITY	SPECIFICITY
Wondfo Dengue NS1 Antigen Test	98.6%	99.1%
Wondfo Dengue IgG/IgM Antibody Test	99.02%	96.8%

TEST PROCEDURES

For Dengue NS1 Antigen Test (15-20 mins.)



For Dengue IgG/IgM Antigen Test (15-20 mins.)



Also available in Dengue IgG/IgM Antibody and Dengue NS1 Antigen Test.
All items are available in 25 test/ kit.