

HCG PREGNANCY URINE / BLOOD TEST

Wondfo Pregnancy Test is a rapid pregnancy, which you can easily carry out by yourself. It detects the presence Human Chronic Gonadotropin (hCG), which appears in urine very early during pregnancy.

Humanitarian Design:





Important:

Do not allow the urine level to exceed the Max (marker line), otherwise the test may not perform correctly. The complete reaction time of 3 minutes is required and your results should appear within 3 minutes. Do not read results after 5 minutes.



General Information:

- Analytical Sensitivity: 25mlU/ml of hCG
- Sensitivity: 99.9%, Specificity: 99.9%
- Specimen: Urine/Serum
- Reading Time: 3 mins. for Urine, 5 mins. for Serum
- Immunochromatography based one-step in vitro diagnostic test
- No cross reactivity: hLH (300mlU/ml), hFSH (300mlU/ml), and hTSH (1,000µlU/ml)
- Storage: 4~30°C

