



AGGLUTINATION VIEWER

Provides a magnified view of the undersurface of test tube contents



FEATURES :

- Designed for use in blood typing, crossmatching, prothrombin time, Kahn and other agglutinations, flocculation, and serodiagnostic tests.
- Illuminates test tube contents through the use of a 5-watt bulb from above and a magnifying mirror from below.
- Easy to use side-mounted ON/OFF switch.
- Flexible neck allows for easy manipulation of light source.
- Input: 120 - 220V 50/60 Hz, 18w; Output: 12VAC 500mA.

SPECIFICATIONS :

Certifications/Compliance	UL Listed
Includes	Bulb, magnifying mirror, power pack, side mounted on/off switch, flexible neck
Wattage	5 Watts
Hertz	50/60Hz
Voltage	120 - 220V