



## CLINICAL CENTRIFUGE DM0408

### FEATURES:

- Speed range of 300-4000rpm.
- 3 programs storage capacity for routine applications.
- Automatic internal diagnosis.
- Maintenance free brushless DC motor.



### ROTORS:

#### A8-15P

Max. Speed: 4000rpm  
Rotor Capacity: 15ml×8  
Rotor Material: High strength plastic  
Default Adaptor: A15P22



#### A12-10P

Rotor Capacity: 15ml×8, 10ml×12  
Rotor Material: High strength plastic

### ADAPTORS:

#### A15P22

15ml adapter for 17×120mm centrifuge tube, 13/16×100mm blood collection tube  
24 pcs/pk



#### A2P17

1.5/2ml Cryogenic Vials/ centrifuge tube adapter  
24 pcs/pk



### SPECIFICATIONS

### DM0408

Speed Range	300-4000rpm : 100 rpm
Max. RCF	1900 x g
Rotor	Angle 8 x 15ml (Default) Angle 10ml x 12
Run Time	1-99 min, HOLD
Motor	Brushless DC Motor
Display	LED
Program	Blood: 3200 rpm, 10 min. Urine: 1800 rpm, 5 min. Fecal: 1300 rpm, 10 min.
Safety Device	Automatic Internal Diagnosis
Acceleration / Deceleration	25s↑ /35s↓
Power	Single-Phase 220V-240V, 50Hz/60Hz
Dimension (W x D x H)	286 mm x 367 mm x 227 mm
Weight	9.2 Kg
Noise	≤60dB(A)



**TruLaboratories**  
Corporation

+632 82-TRULAB (878522)  
+63917.878.2288 • +63925.878.2288  
www.trulaboratories.com