



VORTEX MIXER MX-S

- Touch operation or continuous mode
- Variable speed control from 0 to 3000rpm
- Used for various mixing applications with optional adapters
- Specially designed vacuum suction feet for body stability
- Robust aluminum-cast construction



SPECIFICATIONS	
Voltage	220V, 50/60Hz
Power	60W
Mixing Motion	Orbital
Orbital Diameter	4mm
Motor Type	Shaded Pole Motor
Motor Rating Input	58W
Motor Rating Output	10W
Speed Range	100-3000rpm
Speed Display	Scale
Run Type	Touch Operation / Continuous
Dimensions [WxHxD]	127 x 130 x 160mm
Weight	3.5 Kg.
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	80% RH
Protection Class	IP21

TUBE ADAPTORS	REMARKS
	48 holes test tubes Ø6mm
	15 holes test tubes Ø10mm
	12 holes test tubes, Ø12mm
	For 8 holes test tubes Ø16mm
	For 8 holes test tubes Ø20mm

ACCESSORIES (OPTIONAL)



Standard Top (Default)
for \varnothing30mm tubes and small vessels



Tube Holding Rod
used with tube adapters



Universal Top Plate
Ø100mm



Platform Pad
\varnothing99mm tubes and small vessels



Vacuum Chuck
Made of Rubber



+632 82-TRULAB (878522)
+63917.878.2288 • +63925.878.2288
www.trulaboratories.com