

## CLINICAL CENTRIFUGE DM0408: 12 Placer with 3Pro

DLAB low speed clinical centrifuges are widely used in the medical industry. Ideal for separation of serum, plasma, urine, fecal samples in biochemical & medical laboratories, it is an essential equipment of every clinical laboratory, blood bank, medical college and medical research institute.



### Features

- Speed range of 300-4000rpm.
- 3 programs storage capacity for routine applications.
- Automatic internal diagnosis.
- Maintenance free brushless DC motor.

### Specifications

Speed Range	300-4000rpm, increment: 100rpm
Max. RCF	1900Xg
Speed accuracy	±100rpm
Rotor	Angel 8×15ml (Default), Angel 10ml×12
Run Time	1-99min, HOLD
Motor	Brushless DC motor
Display	LED
Program	Blood: 3200rpm, 10min. Urine: 1800rpm, 5min. Faecal: 1300rpm, 10min.
Safety Devices	Automatic internal diagnosis
Acceleration/Braking time	25s↑35s↓
Power	Single phase, 100V-240V, 50Hz/60Hz
Dimension [D×W×H]	286×367×227mm
Weight	9.2kg
Noise	≤60dB(A)

