

ROTATOR & DIGITAL ROTATOR

DSR-Series



DSR-2100P



DSR-2800P



DSR-2100D



DSR-2800D

Features:

- * Wide used in VDRL, PRP, and EIA test.
- * Stable and quiet operation.
- * Simple and digital type for your options.

- * With variable speed, timer control, and continuous operation.
- * Quality certification: ISO 9001, ISO 13485

Specifications:

Model	DSR-2100P	DSR-2800P	DSR-2100D	DSR-2800D
Speed	20~240rpm			
Timer	0~30 min. mechanical timer		99 hours and 59 minutes digital timer	
Display	No display		Digital display for present speed value and remaining time	
Platform dimension	230x230(mm)	280x280(mm)	230x230(mm)	280x280(mm)
Orbit diameter	20mm			
Load capacity	2.0kgs at the speed from 20~150rpm; 1.0kgs at the speed from 151~240rpm			
Accessory	1 pc of spring			
Weight (kg) (N.W./G.W.)	3.7 / 4.5	4.0 / 4.8	3.9 / 4.7	4.2 / 5.0
Overall dimension (mm)(WxDxH)	260x295x140	280x295x140	260x295x140	280x295x140

Power supply: AC 110V 60Hz (220V 50/60Hz available)

ORBITAL SHAKER

OS-Series



OS-500 with flat platform



OS-500 with clamps

Features:

- * Designed for continuous-operation orbital shaking, day and night.
- * Suitable for carrying heavy duty.
- * Maintenance-free induction brushless motor.
- * Differential flask clamps available.
- * All stainless steel platform and clamps.
- * Quality certification: ISO 9001, ISO 13485

Specifications:

Model	OS-340	OS-500	OS-600	OS-1060
Speed	10~250rpm at 60Hz; 10~240rpm at 50Hz			
Platform dimension	W340xD340(mm)	W500xD500(mm)	W600xD600(mm)	W1060xD600(mm)
Orbit diameter	25mm			
Optional Accessory	50ml, 125ml, 250ml, 500ml, 1L, 2L flask clamps			
Overall dimension (mm)(WxDxH)	380x410x195	560x600x225	660x580x225	1120x700x225

Power supply: AC 110V 60Hz (220V 50/60Hz available)



DIGISYSTEM LABORATORY INSTRUMENTS INC.

NO. 1, LANE 393, FU DE 1 ROAD, XIZHI DIST., NEW TAIPEI CITY, TAIWAN
 TEL: 886-2-26941738 . 26942339 . 26945547 FAX: 886-2-26941695
 E-MAIL: digilab@ms23.hinet.net INTERNET: www.digisystemlab.com
 We reserve the rights to change the specifications without notice.

V.201301