

ANALOG CENTRIFUGE

FOR CLINICAL AND RESEARCH APPLICATIONS

FEATURES:

- Operation can be timed from 0 seconds to 30 minutes or continuous running
- Gentle breaking low speeds with efficient separation
- Quick spins are possible by pressing and holding PULSE key
- High strength plastic rotor and excellent balance technology for quiet and stable operation



Max. Speed [rpm]	0-4500 rpm, Analog Control
Max. RCF [xg]	2490xg
Speed Accuracy [rpm]	±20
Rotor Capacity	12x10ml/7ml/5ml; 12x15ml
Run Time	0 seconds to 30 minutes Timer or Continuous Operation
Driving Motor	DC motor Brush
Weight	14 pounds
Power	110V-240V, 50Hz/60Hz, 3A
Dimensions	(L) 11.0 x (W) 10 x (H) 5.5 inches
Weight	6 kg

