

TruLaboratories COVID-19 IgG/IgM Rapid Tests

Wondfo and Innovita: A Quick Comparison Guide

Wondfo®

One box contains:

- 20 test kits
- 20 blood pipettes
- 1 buffer solution
- 1 package insert



Total IgG/IgM Antibody Detection

Assay Performance

Sensitivity	86.43%
Specificity	99.57%

INNOVITA™

One box contains:

- 40 test kits
- 1 buffer solution
- 1 package insert



Separate indicators for IgG and IgM

Assay Performance

Sensitivity	91.39%
Specificity	98.45%

Recommended Sample Collection Method



Using the venous whole blood, serum, or plasma from the tube

- Negative results do not completely eliminate the COVID-19 infection, particularly in those who have been in contact with the virus.
- Undergo confirmatory testing with RT-PCR to rule out infection in individuals.
- Results from antibody testing should not be the sole basis to diagnose and confirm an infection.

Ordering Information

Rapid Test Kits	Quantity
Wondfo COVID-19	20 tests
Innovita COVID-19	40 tests

Other Essential Products for Venipuncture	Quantity
Phlebotomy Arm Support	1 unit
Multisample Needle	100 pcs/box
Needle Adaptor	1 pc
Tourniquet	1 pc
TruLab EDTA Blood Tube (Whole Blood)	100 pcs/box
TruLab Clot Activator Blood Tube (Serum)	100 pcs/box

Test Preparation and Safety Disposal	Quantity
Vitrex Alcohol Prep Pads	100 pcs/box
Sharps Container	1 unit

Other Essential Products for Finger Prick	Quantity
Lancing Device	1 unit
Blood Lancet	100 pcs/box
Vitrex Safety Blood Lancets (No Lancing Device Needed)	100 pcs/box



14 Detroit St. Cubao, Quezon City
 8451.1964 | 8451.1974 | 8451.1921
 0917 878.2288 | 0925 878 2288
www.trulaboratories.com