



STANDARD d-BLOCK Incubator is an auxiliary device providing a constant temperature during the test. This product is designed for IVD products that require thermal incubation.

SPECIFICATIONS	
Dimension	220*184*73mm
Initial time	15 minutes
Set Temperature Range	35~40°C / 95~104°F
Accuracy of Temperature	+/- 1°C
Environment Condition	Temperature: 10° C ~ 30° C (50° F to 86° F) / Humidity: 20% ~ 80% Non condensing
Storage Condition	Temperature: 0° C ~ 70° C (32° F to 125° F) / Humidity: $10 \sim 90\%$
Equipment Control	4 buttons
Equipment Measurement Unit	°C, °F
Equipment Display Type	LCD (Customized)
Weight	1.9 Kg
Equipment Ratings	12V (DC), 5A



