

TruScreen

PREGNANCY AND DRUG TEST CONTROLS



PREGNANCY LIQUID CONTROLS

TruScreen Human Chorionic Gonadotropin (hCG) Positive and Negative controls are prepared from human based urine. These hCG controls can be used to monitor the performance of qualitative immunochromatographic hCG detection devices.

- TruScreen hCG Liquid Control Urine, Positive, 5 mL vial
- TruScreen hCG Liquid Control Urine, Negative, 5 mL vial


DRUG TEST LIQUID CONTROL URINE

TruScreen Drug Test Liquid control is an In Vitro Diagnostic (IVD) device, that is intended for use as quality control urine to monitor the precision of laboratory urine toxicology testing procedures for the analytes listed below

TEST	SAMHSA MANDATED	Negative	2X
THC	Delta-9-THC-COOH	0	100
COC	Benzoyllecgonine	0	600
PCP	Phencyclidine	0	50
OPI/MOR	Morphine (Low Opiate)	0	600
	Morphine (High Opiate)	0	4000
AMP	d-Amphetamine	0	2000
MET	d-Methamphetamine	0	2000
NON-MANDATED			
BARB	Secobarbital	0	600
BZO	Oxazepam	0	600
MTD	Methadone	0	600
MQL	Methaqualone	0	600
PPX	Propoxyphene	0	600
TCA	Nortriptyline	0	2000
MDMA	MDMA*	0	1000
OXY	Oxycodone	0	200
BUP	Buprenorphine	0	20



- TruScreen with MDMA, OXY & BUP, 2X Cutoff, 5 mL vial
- TruScreen Negative Control Urine, 5 mL vial



Made In U S A

EMERGO EUROPE

Molenstraat 15

2513 BH, The Hague

The Netherlands



14 Detroit Street, Cubao, Quezon City 1111

Phone: +63 2 451 1964 ; 451 1974 | Fax: +63 2 451 1921

Email: sales@trulaboratories.com

www.trulaboratories.com