




SARS-CoV-2 Antibody Test (Lateral Method)

- As an auxiliary and supplementary method for the further diagnosis to the patient with negative result in PCR.
- Easy to use, instant result in 15 minutes

One box contains:

- 20 test kits
- 20 blood pipettes
- 1 buffer solution
- 1 package insert



Antibodies can be detected as early as one week after infection, with slightly difference between individuals.

1. Negative result indicates no infection or infection at a very early stage (within first few days)
2. Positive results indicates an infection with coronavirus.

Total IgG/IgM Antibody Detection



Assay Performance


Sensitivity	86.43%
Specificity	99.57%

Operation procedure




Fingerstick
Whole Blood

OR




Venous Whole
Blood/Serum/Plasma

01.



02.



03.



15 minutes

04.

