

On•Call® Vivid

Premium Features At An Affordable Price



- Fast, Reliable Results
- User - Friendly for Easy Daily Testing
- Advanced Features for Improving Patient Outcomes

Feature	Specification
Measurement Range	20– 600 mg/dL
Result Calibration	Plasma-equivalent
Sample	Fresh capillary whole blood
Minimum Sample Size	0.8 µL
Test Time	5 seconds
On/Off Source	Two (2) CR 2032 3.0V coin cell battery
Battery Life	1,000 tests
Glucose Units of Measure	The meter is pre-set to milligrams per deciliter (mg/dL)
Memory	Up to 500 records with time and date
Automatic shutoff	2 minutes after last action
Meter Size	89.6 mm x 58 mm x 21.7 mm
Display Size	37 mm x 33.5 mm
Operating Temperature	41 – 113 °F
Operating Relative Humidity	10 – 90%
Hematocrit Range	20 – 70%
Data Communications	USB Cable



User - Friendly for Easy Daily Testing

- Backlight and strip port light allow user to test in the dark.
- Strip ejector allows for fast and easy disposal to avoid contact with used strips.
- Lancet ejector designed for safe and convenient lancet release.

Ordering Information:

TRULABORATORIES CORPORATION

14 Detroit Street, Cubao, Quezon City

Tel. No.: +63 2 451 1964, 451 1974

Fax No.: +63 2 451 1921

Mobile : +63 917 878 2288, 925 878 2288

E-mail : sales@trulaboratories.com



ACON Laboratories, Inc., San Diego, CA 92121, U.S.A.